

HSP ?2 Polyclonal Antibody

Cat No: HR1AP3690

For research use only

Overview

Product Name	HSP ?2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	HSPB2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HSPB2; Heat shock protein beta-2; HspB2; DMPK-binding protein; MKBP
Background	heat shock protein family B (small) member 2(HSPB2) Homo sapiens The protein encoded by this gene belongs to the superfamily of small heat-shock proteins containing a conservative alpha-crystallin domain at the C-terminal part of the molecule. The protein is expressed preferentially in the heart and skeletal muscle. This protein regulates Myotonic Dystrophy Protein Kinase, which plays an important role in maintenance of muscle structure and function. [provided by RefSeq, Dec 2012],