

GALK1 Polyclonal Antibody

Cat No: HR1AP3507

For research use only

Overview

Product Name	GALK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	GALK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GALK1; GALK; Galactokinase; Galactose kinase
Background	galactokinase 1(GALK1) Homo sapiens Galactokinase is a major enzyme for the metabolism of galactose and its deficiency causes congenital cataracts during infancy and presenile cataracts in the adult population. [provided by RefSeq, Jul 2008],