

dTMP Kinase Polyclonal Antibody

Cat No: HR1AP3316

For research use only

Overview

Product Name	dTMP Kinase Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	DTYMK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DTYMK; CDC8; TMPK; TYMK; Thymidylate kinase; dTMP kinase
Background	catalytic activity:ATP + dTMP = ADP + dTDP.,function:Catalyzes the conversion of dTMP to dTDP.,pathway:Pyrimidine metabolism; dTTP biosynthesis.,similarity:Belongs to the thymidylate kinase.,