

## CLN5 Polyclonal Antibody

Cat No: HR1AP5793

For research use only

### Overview

Product Name	CLN5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	CLN5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLN5; Ceroid-lipofuscinosis neuronal protein 5; Protein CLN5
Background	ceroid-lipofuscinosis, neuronal 5(CLN5) Homo sapiens This gene is one of eight which have been associated with neuronal ceroid lipofuscinoses (NCL). Also referred to as Batten disease, NCL comprises a class of autosomal recessive, neurodegenerative disorders affecting children. The genes responsible likely encode proteins involved in the degradation of post-translationally modified proteins in lysosomes. The primary defect in NCL disorders is thought to be associated with lysosomal storage function.[provided by RefSeq, Oct 2008].