

CEP78 Polyclonal Antibody

Cat No: HR1AP3064

For research use only

Overview

Product Name	CEP78 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	CEP78
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEP78; C9orf81; Centrosomal protein of 78 kDa; Cep78
Background	similarity:Belongs to the CEP78 family.,similarity:Contains 3 LRR (leucine-rich) repeats.,