

CEP70 Polyclonal Antibody

Cat No: HR1AP3063

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | CEP70 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Monkey |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 70kDa |
| GeneID?Human? | CEP70 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | CEP70; BITE; Centrosomal protein of 70 kDa; Cep70; p10-binding protein |
| Background | similarity:Contains 1 TPR repeat., |