

CDKAL1 Polyclonal Antibody

Cat No: HR1AP7792

For research use only

Overview

Product Name	CDKAL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	CDKAL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDKAL1; Threonylcarbamoyladenine tRNA methylthiotransferase; CDK5 regulatory subunit-associated protein 1-like 1; tRNA-t(6)A37 methylthiotransferase
Background	CDK5 regulatory subunit associated protein 1 like 1(CDKAL1) Homo sapiens The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes. [provided by RefSeq, May 2010],