

# CEP152 Polyclonal Antibody

Cat No: HR1AP6396

For research use only

## Overview

Product Name	CEP152 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	CEP152
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEP152; KIAA0912; Centrosomal protein of 152 kDa; Cep152
Background	centrosomal protein 152(CEP152) Homo sapiens This gene encodes a protein that is thought to be involved with centrosome function. Mutations in this gene have been associated with primary microcephaly (MCPH4). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2010],