

C1qL2 Polyclonal Antibody

Cat No: HR1AP2920

For research use only

Overview

Product Name	C1qL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32-36kDa
GeneID?Human?	C1QL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C1QL2; Complement C1q-like protein 2
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,