

Arylsulfatase I Polyclonal Antibody

Cat No: HR1AP2822

For research use only

Overview

Product Name	Arylsulfatase I Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	ARSI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARSI; Arylsulfatase I; ASI
Background	arylsulfatase family member I(ARSI) Homo sapiens This gene encodes a protein that belongs to a large family of sulfatases that hydrolyze sulfate esters and sulfamates. Members of this family play a role in several cellular processes, including hormone synthesis, cell signaling in development and degradation of macromolecules. The protein encoded by this gene is thought to be secreted, and to function in extracellular space. [provided by RefSeq, Jul 2016],