

ARMC6 Polyclonal Antibody

Cat No: HR1AP2814

For research use only

Overview

Product Name	ARMC6 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	ARMC6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARMC6; Armadillo repeat-containing protein 6
Background	armadillo repeat containing 6(ARMC6) Homo sapiens The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],