



Annexin X Polyclonal Antibody

Cat No: HR1AP2770

For research use only

Overview

Product Name	Annexin X Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	ANXA10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14
Background	annexin A10(ANXA10) Homo sapiens This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],