

HSP90? Rabbit Polyclonal Antibody

Cat No: HR1AP9671

For research use only

Overview

Product Name	HSP90? Rabbit Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	HSP90AA1 HSP90A HSPC1 HSPCA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)
Background	heat shock protein 90 alpha family class A member 1(HSP90AA1) Homo sapiens The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],