

# ZBTB45 Polyclonal Antibody

Cat No: HR1AP9289

For research use only

## Overview

Product Name	ZBTB45 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	70kDa
GeneID?Human?	ZBTB45
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Zinc finger and BTB domain-containing protein 45; Zinc finger protein 499
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,