

## PAR-6? Polyclonal Antibody

Cat No: HR1AP8876

For research use only

### Overview

Product Name	PAR-6? Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	PARD6B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PARD6B; PAR6B; Partitioning defective 6 homolog beta; PAR-6 beta; PAR-6B
Background	<p>par-6 family cell polarity regulator beta(PARD6B) Homo sapiens This gene is a member of the PAR6 family and encodes a protein with a PSD95/Discs-large/ZO1 (PDZ) domain, an OPR domain and a semi-Cdc42/Rac interactive binding (CRIB) domain. This cytoplasmic protein is involved in asymmetrical cell division and cell polarization processes as a member of a multi-protein complex. [provided by RefSeq, Jul 2008],</p>