

NDRG2 Polyclonal Antibody

Cat No: HR1AP8007

For research use only

Overview

Product Name	NDRG2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	NDRG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDRG2; KIAA1248; SYLD; Protein NDRG2; N-myc downstream-regulated gene 2 protein; Protein Syld709613
Background	NDRG family member 2(NDRG2) Homo sapiens This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2008],