



eIF1B Polyclonal Antibody

Cat No: HR1AP5566

For research use only

Overview

Product Name	eIF1B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	EIF1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EIF1B; Eukaryotic translation initiation factor 1b; eIF1b; Protein translation factor SUI1 homolog GC20
Background	function:Probably involved in translation.,similarity:Belongs to the SUI1 family.,