

ZAR1 Polyclonal Antibody

Cat No: HR1AP4834

For research use only

Overview

Product Name	ZAR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	ZAR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZAR1; Zygote arrest protein 1; Oocyte-specific maternal effect factor
Background	zygote arrest 1(ZAR1) Homo sapiens This maternal effect gene is oocyte-specific and encodes a protein that is thought to function in the initiation of embryogenesis. A similar protein in mouse is required for female fertility. [provided by RefSeq, Jul 2013],