

# Neurensin-1 Polyclonal Antibody

Cat No: HR1AP4037

For research use only

## Overview

Product Name	Neurensin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	NRSN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NRSN1; VMP; Neurensin-1; Neuro-p24; Vesicular membrane protein of 24 kDa; Vesicular membrane protein p24
Background	function:May play an important role in neural organelle transport, and in transduction of nerve signals or in nerve growth. May play a role in neurite extension.,similarity:Belongs to the VMP family.,subcellular location:Localizes mainly to neurites.,tissue specificity:Expressed in brain. Not detectable in other tissues tested.,