

MOV10L1 Polyclonal Antibody

Cat No: HR1AP3925

For research use only

Overview

Product Name	MOV10L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	MOV10L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MOV10L1; Putative helicase Mov10l1; Moloney leukemia virus 10-like protein 1; MOV10-like protein 1
Background	Mov10 RISC complex RNA helicase like 1(MOV10L1) Homo sapiens This gene is similar to a mouse gene that encodes a putative RNA helicase and shows testis-specific expression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009].