

Ataxin-2L Polyclonal Antibody

Cat No: HR1AP5739

For research use only

Overview

Product Name	Ataxin-2L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	113kDa
GeneID?Human?	ATXN2L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATXN2L; A2D; A2LG; A2LP; A2RP; Ataxin-2-like protein; Ataxin-2 domain protein; Ataxin-2-related protein
Background	ataxin 2 like(ATXN2L) Homo sapiens This gene encodes an ataxin type 2 related protein of unknown function. This protein is a member of the spinocerebellar ataxia (SCAs) family, which is associated with a complex group of neurodegenerative disorders. Several alternatively spliced transcripts encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],