

# Erbin (phospho Tyr1104) Polyclonal Antibody

Cat No: HR1AP2418

For research use only

## Overview

Product Name	Erbin (phospho Tyr1104) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	172kDa
GeneID?Human?	ERBB2IP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ERBB2IP; ERBIN; KIAA1225; LAP2; Protein LAP2; Densin-180-like protein; Erb2-interacting protein; Erbin
Background	<p>erbb2 interacting protein(ERBIN) Homo sapiens This gene is a member of the leucine-rich repeat and PDZ domain (LAP) family. The encoded protein contains 17 leucine-rich repeats and one PDZ domain. It binds to the unphosphorylated form of the ERBB2 protein and regulates ERBB2 function and localization. It has also been shown to affect the Ras signaling pathway by disrupting Ras-Raf interaction. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011],</p>