

# Dok-2 (phospho Tyr299) Polyclonal Antibody

Cat No: HR1AP2410

For research use only

## Overview

Product Name	Dok-2 (phospho Tyr299) Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	DOK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DOK2; Docking protein 2; Downstream of tyrosine kinase 2; p56(dok-2)
Background	docking protein 2(DOK2) Homo sapiens The protein encoded by this gene is constitutively tyrosine phosphorylated in hematopoietic progenitors isolated from chronic myelogenous leukemia (CML) patients in the chronic phase. It may be a critical substrate for p210(bcr/abl), a chimeric protein whose presence is associated with CML. This encoded protein binds p120 (RasGAP) from CML cells. [provided by RefSeq, Jul 2008],