

## Septin 3 Polyclonal Antibody

Cat No: HR1AP4536

For research use only

### Overview

Product Name	Septin 3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	Sep-03
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SEPT3; SEP3; Neuronal-specific septin-3
Background	<p>septin 3(SEPT3) Homo sapiens This gene belongs to the septin family of GTPases. Members of this family are required for cytokinesis. Expression is upregulated by retinoic acid in a human teratocarcinoma cell line. The specific function of this gene has not been determined. Alternative splicing of this gene results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],</p>