

4.1G Polyclonal Antibody

Cat No: HR1AP2660

For research use only

Overview

Product Name	4.1G Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	112kDa
GeneID?Human?	EPB41L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPB41L2; Band 4.1-like protein 2; Generally expressed protein 4.1; 4.1G
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 FERM domain.,subunit:The CTD domain interacts with FKBP2 (By similarity). Interacts with FCGR1A.,tissue specificity:Widely expressed.,