

Flotillin-2 Polyclonal Antibody

Cat No: HR1AP9397

For research use only

Overview

Product Name	Flotillin-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49,58kDa
GeneID?Human?	FLOT2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Flotillin-2 (Epidermal surface antigen) (ESA) (Membrane component chromosome 17 surface marker 1)
Background	flotillin 2(FLOT2) Homo sapiens Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. This gene encodes a caveolae-associated, integral membrane protein, which is thought to function in neuronal signaling. [provided by RefSeq, Jul 2008],