

Stathmin-3 Polyclonal Antibody

Cat No: HR1AP4631

For research use only

Overview

Product Name	Stathmin-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	STMN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	STMN3; SCLIP; Stathmin-3; SCG10-like protein
Background	<p>stathmin 3(STMN3) Homo sapiens This gene encodes a protein which is a member of the stathmin protein family. Members of this protein family form a complex with tubulins at a ratio of 2 tubulins for each stathmin protein. Microtubules require the ordered assembly of alpha- and beta-tubulins, and formation of a complex with stathmin disrupts microtubule formation and function. A pseudogene of this gene is located on chromosome 22. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013],</p>