

MSPD3 Polyclonal Antibody

Cat No: HR1AP13099

For research use only

Overview

Product Name	MSPD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	MOSPD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>motile sperm domain containing 3(MOSPD3) Homo sapiens This gene encodes a multi-pass membrane protein with a major sperm protein (MSP) domain. The deletion of a similar mouse gene is associated with defective cardiac development and neonatal lethality. Alternate transcriptional splice variants, encoding different isoforms, have been described. [provided by RefSeq, Jul 2008],</p>