

## COPT1 Polyclonal Antibody

Cat No: HR1AP13094

For research use only

### Overview

Product Name	COPT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	SLC31A1 COPT1 CTR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 31 member 1(SLC31A1) Homo sapiens The protein encoded by this gene is a high-affinity copper transporter found in the cell membrane. The encoded protein functions as a homotrimer to effect the uptake of dietary copper. [provided by RefSeq, Aug 2011],