

MLP3C Polyclonal Antibody

Cat No: HR1AP13088

For research use only

Overview

Product Name	MLP3C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	MAP1LC3C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	microtubule associated protein 1 light chain 3 gamma(MAP1LC3C) Homo sapiens Autophagy is a highly regulated bulk degradation process that plays an important role in cellular maintenance and development. MAP1LC3C is an ortholog of the yeast autophagosome protein Atg8 (He et al., 2003 [PubMed 12740394]).[supplied by OMIM, Nov 2010],