

OSR1 Polyclonal Antibody

Cat No: HR1AP13082

For research use only

Overview

Product Name	OSR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	OSR1 ODD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the Odd C2H2-type zinc-finger protein family.,similarity:Contains 3 C2H2-type zinc fingers.,tissue specificity:Expressed in adult colon, small intestine, prostate, testis, and fetal lung.,