

DERL2 Polyclonal Antibody

Cat No: HR1AP13079

For research use only

Overview

Product Name	DERL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	DERL2 DER2 FLANA CGI-101 SBB153
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>derlin 2(DERL2) Homo sapiens Proteins that are unfolded or misfolded in the endoplasmic reticulum (ER) must be refolded or degraded to maintain the homeostasis of the ER. DERL2 is involved in the degradation of misfolded glycoproteins in the ER (Oda et al., 2006 [PubMed 16449189]).[supplied by OMIM, Mar 2008],</p>