

RBM10 Polyclonal Antibody

Cat No: HR1AP13077

For research use only

Overview

Product Name	RBM10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	RBM10 DXS8237E GPATC9 GPATCH9 KIAA0122
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RNA binding motif protein 10(RBM10) Homo sapiens This gene encodes a nuclear protein that belongs to a family proteins that contain an RNA-binding motif. The encoded protein associates with hnRNP proteins and may be involved in regulating alternative splicing. Defects in this gene are the cause of the X-linked recessive disorder, TARP syndrome. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Mar 2011],</p>