

# SUN5 Polyclonal Antibody

Cat No: HR1AP13072

For research use only

## Overview

Product Name	SUN5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	SUN5 SPAG4L TSARG4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 Unc84 (SUN) domain.,tissue specificity:Widely expressed.,