

## ZEB1 Polyclonal Antibody

Cat No: HR1AP13068

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | ZEB1 Polyclonal Antibody  |
| Source                | Rabbit  |
| Applications          | WB,ELISA  |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 123kDa  |
| GeneID?Human?         | ZEB1 AREB6 TCF8   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     |   |
| Background            | zinc finger E-box binding homeobox 1(ZEB1) Homo sapiens This gene encodes a zinc finger transcription factor. The encoded protein likely plays a role in transcriptional repression of interleukin 2. Mutations in this gene have been associated with posterior polymorphous corneal dystrophy-3 and late-onset Fuchs endothelial corneal dystrophy. Alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2010], |