

## DLL4 Polyclonal Antibody

Cat No: HR1AP13067

For research use only

### Overview

Product Name	DLL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	DLL4 UNQ1895/PRO4341
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	delta like canonical Notch ligand 4(DLL4) Homo sapiens This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. [provided by RefSeq, Jul 2008],