

## BIRC2 Polyclonal Antibody

Cat No: HR1AP13066

For research use only

### Overview

Product Name	BIRC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	BIRC2 API1 IAP2 MIHB RNF48
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>baculoviral IAP repeat containing 2(BIRC2) Homo sapiens The protein encoded by this gene is a member of a family of proteins that inhibits apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. This encoded protein inhibits apoptosis induced by serum deprivation and menadione, a potent inducer of free radicals. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],</p>