

VGLU2 Polyclonal Antibody

Cat No: HR1AP13065

For research use only

Overview

Product Name	VGLU2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	SLC17A6 DNPI VGLUT2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	developmental stage:Expressed in fetal brain.,function:Mediates the uptake of glutamate into synaptic vesicles at presynaptic nerve terminals of excitatory neural cells. May also mediate the transport of inorganic phosphate.,similarity:Belongs to the major facilitator superfamily. Sodium/anion cotransporter family. VGLUT subfamily.,tissue specificity:Predominantly expressed in adult brain. Expressed in amygdala, caudate nucleus, cerebral cortex, frontal lobe, hippocampus, medulla, occipital lobe, puteman, spinal cord, substantia nigra, subthalamic nucleus, temporal lobe and thalamus.,