

MIEN1 Polyclonal Antibody

Cat No: HR1AP13061

For research use only

Overview

Product Name	MIEN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	MIEN1 C17orf37 RDX12 XTP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in a redox-related process .,similarity:Belongs to the SeiWTH family.,