

VAMP7 Polyclonal Antibody

Cat No: HR1AP13055

For research use only

Overview

Product Name	VAMP7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	VAMP7 SYBL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	vesicle associated membrane protein 7(VAMP7) Homo sapiens This gene encodes a transmembrane protein that is a member of the soluble N-ethylmaleimide-sensitive factor attachment protein receptor (SNARE) family. The encoded protein localizes to late endosomes and lysosomes and is involved in the fusion of transport vesicles to their target membranes. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jun 2010],