

RFA4 Polyclonal Antibody

Cat No: HR1AP13047

For research use only

Overview

Product Name	RFA4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	RPA4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>replication protein A4(RPA4) Homo sapiens This gene encodes a single-stranded DNA-binding protein that is the 30-kDa subunit of the replication protein A complex. Replication protein A is an essential factor for DNA double-strand break repair and cell cycle checkpoint activation. The encoded protein localizes to DNA repair foci and may be involved in the cellular DNA damage response. This protein may also play a role in inhibiting viral replication.[provided by RefSeq, Apr 2010],</p>