

PLXA4 Polyclonal Antibody

Cat No: HR1AP13038

For research use only

Overview

Product Name	PLXA4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	208kDa
GeneID?Human?	PLXNA4 KIAA1550 PLXNA4A PLXNA4B UNQ2820/PRO34003
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Forms a receptor complex with neuropilin-1 and can propagate semaphorin-3A elicited inhibitory signals into cells and neurons.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the plexin family.,similarity:Contains 1 Sema domain.,similarity:Contains 3 PSI domains.,similarity:Contains 4 IPT/TIG domains.,