

PCNT Polyclonal Antibody

Cat No: HR1AP13035

For research use only

Overview

Product Name	PCNT Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	366kDa
GeneID?Human?	PCNT KIAA0402 PCNT2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>pericentrin(PCNT) Homo sapiens The protein encoded by this gene binds to calmodulin and is expressed in the centrosome. It is an integral component of the pericentriolar material (PCM). The protein contains a series of coiled-coil domains and a highly conserved PCM targeting motif called the PACT domain near its C-terminus. The protein interacts with the microtubule nucleation component gamma-tubulin and is likely important to normal functioning of the centrosomes, cytoskeleton, and cell-cycle progression. Mutations in this gene cause Seckel syndrome-4 and microcephalic osteodysplastic primordial dwarfism type II. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015],</p>