

## ACPM Polyclonal Antibody

Cat No: HR1AP13030

For research use only

### Overview

Product Name	ACPM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	NDUFAB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Carrier of the growing fatty acid chain in fatty acid biosynthesis in mitochondria. Accessory and non-catalytic subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), which functions in the transfer of electrons from NADH to the respiratory chain.,similarity:Contains 1 acyl carrier domain.,subunit:Mammalian complex I is composed of 45 different subunits.,