

ATOH7 Polyclonal Antibody

Cat No: HR1AP13028

For research use only

Overview

Product Name	ATOH7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	ATOH7 ATH5 BHLHA13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>atonal bHLH transcription factor 7(ATOH7) Homo sapiens This intronless gene encodes a member of the basic helix-loop-helix family of transcription factors, with similarity to Drosophila atonal gene that controls photoreceptor development. Studies in mice suggest that this gene plays a central role in retinal ganglion cell and optic nerve formation. Mutations in this gene are associated with nonsyndromic congenital retinal nonattachment. [provided by RefSeq, Dec 2011],</p>