

STAR3 Polyclonal Antibody

Cat No: HR1AP13029

For research use only

Overview

Product Name	STAR3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	STARD3 CAB1 MLN64
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>STAR related lipid transfer domain containing 3(STARD3) Homo sapiens This gene encodes a member of a subfamily of lipid trafficking proteins that are characterized by a C-terminal steroidogenic acute regulatory domain and an N-terminal metastatic lymph node 64 domain. The encoded protein localizes to the membranes of late endosomes and may be involved in exporting cholesterol. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009],</p>