

HEY1 Polyclonal Antibody

Cat No: HR1AP13019

For research use only

Overview

Product Name	HEY1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	HEY1 BHLHB31 CHF2 HERP2 HESR1 HRT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>hes related family bHLH transcription factor with YRPW motif 1(HEY1) Homo sapiens This gene encodes a nuclear protein belonging to the hairy and enhancer of split-related (HESR) family of basic helix-loop-helix (bHLH)-type transcriptional repressors. Expression of this gene is induced by the Notch and c-Jun signal transduction pathways. Two similar and redundant genes in mouse are required for embryonic cardiovascular development, and are also implicated in neurogenesis and somitogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008],</p>