



# HIC2 Polyclonal Antibody

Cat No: HR1AP13020

For research use only

## Overview

Product Name	HIC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	HIC2 HRG22 KIAA1020 ZBTB30
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family. Hic subfamily.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Interacts with CtBP.,tissue specificity:Highest levels in cerebellum.,