

S29A2 Polyclonal Antibody

Cat No: HR1AP13011

For research use only

Overview

Product Name	S29A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	SLC29A2 DER12 ENT2 HNP36
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 29 member 2(SLC29A2) Homo sapiens The uptake of nucleosides by transporters, such as SLC29A2, is essential for nucleotide synthesis by salvage pathways in cells that lack de novo biosynthetic pathways. Nucleoside transport also plays a key role in the regulation of many physiologic processes through its effect on adenosine concentration at the cell surface (Griffiths et al., 1997 [PubMed 9396714]).[supplied by OMIM, Nov 2008],</p>