

MARE3 Polyclonal Antibody

Cat No: HR1AP13008

For research use only

Overview

Product Name	MARE3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	MAPRE3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	microtubule associated protein RP/EB family member 3(MAPRE3) Homo sapiens The protein encoded by this gene is a member of the RP/EB family of genes. The protein localizes to the cytoplasmic microtubule network and binds APCL, a homolog of the adenomatous polyposis coli tumor suppressor gene. [provided by RefSeq, Jul 2008],